

CAREMED ALRICK ENDLESS SERIES

SKU: EN7KS-0000

PRODUCT DESCRIPTION

The EN7 series bed is a variant of the EN9000 manufactured overseas. It is the first ever universal bed that suits low care, high care and bariatric care, plus the innovative three-in-one bed. The EN7 uses Linak motors and has a low height of 120mm and has no towers!

SIZE OPTIONS

Available in single and king single models.

BENEFITS

It is the first ever universal bed that suits low care, high care and bariatric care, plus the innovative three-in-one bed. With the extremely powerful Linak motors, this bed truly has no towers! Even with so many features packed into the one bed, controlling the bed is simpler than ever. Tailor this bed to all your needs for every moment of the day and night. No need to stress about the things that were once made hard with bedding, allowing you to get back to what matters most to you. Here are just a few of the features...

FEATURES

Feature	Details
Electric Hi/Low	Auto regression backrest
Electric knee bend	75mm Twin wheel castors
Electric forward & reverse Trendelenburg	Standard head and foot boards (Finl Blonde laminate)

Item	Details

Item	Details
Overall dimensions (mm)	2080 × 920 (Single) 2080 × 1090 (King single)
Sleeping surface dimensions (mm)	2000 × 890 (Single) 2000 × 1070 (King single)
SWL	200kg / 300kg
Net bed weight	110kg
Low position	120mm
High position	725mm
Underbed clearance	Clearview

Function	Degrees
Backrest	62°
Knee	35°
Trendelenburg	15°
Reverse Trendelenburg	15°

Component	Warranty
Frame	10 years
Electrics	5 years
Handset & Castors	1 year

Options

Underbed light
Battery backup
Dual braking
End board options

Brochure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au