

ASPIRE CLASSIC BEDSIDE COMMODE CHAIR

SKU: BEC04500-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

A neat yet discrete design that allows for toileting in any room. Our seat cushions are made with quality materials designed for maximum comfort. This means you can feel happy knowing this is there when you need it yet it can still be used a comfortable resting place. Rest assured our frames are designed to withstand large amounts of weight and force making for a very versatile unit.

FEATURES

- Raised Padded Armrests
- Easy Clean Vinyl Surfaces
- Height Adjustable
- Removable, Padded Seat
- Non-Slip Rubber Tips
- Seat Height: 475 - 615 mm
- Weight Capacity: 125kg
- More information in the 'Specifications' tab.

SPECIFICATIONS

Safe Working Load	125 kg
Total Height	900 - 1040 mm
Total Length	600 mm
Total Width	560 mm

Backrest Height	450 mm
Seat Depth	440 mm
Seat To Floor Height	475 - 615 mm
Seat Width	450 mm
Unit Weight	5.5 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au