

CONFIGURA ADVANCE MANUAL CARE CHAIR

SKU: CR5438

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Configura Advance Manual Care Chair, a highly configurable specialist chair designed to meet complex and challenging seating needs without breaking the bank. This innovative and cost-effective seating solution combines state-of-the-art features and standards that set it apart from other chairs in its class.

FEATURES

- **37 Degrees of Tilt in Space & Full-Length Compensating Legrest:** Experience unparalleled comfort and flexibility with a wide range of tilt options, along with a compensating legrest.
- **Lockable, Reduced Shear Backrest Recline:** Maximize pelvic stability with a backrest recline feature that minimizes shear forces, giving you the support you need.
- **Tool-less Adjustments with ClickLock Adjustment System:** Say goodbye to cumbersome tools and enjoy easy and hassle-free adjustments with the innovative ClickLock system.
- **Upgradeable Backrest Options:** Tailor your chair to your specific needs with customizable backrest options.
- **Removable, Angle, and Height-Adjustable Footplate as Standard:** For personalized comfort and fit, this chair comes with a footplate that can be adjusted in height and angle or entirely removed.
- **Visco Integrated Self-Supply Cushion as Standard:** Measuring 17.5 x 17.5 x 3 inches, the self-supply cushion ensures additional comfort and relaxation.
- **Drop-Down Armrest Option:** Facilitating easier transfers and sling insertions, the drop-down armrest is a functional feature that adds to the chair's versatility.

- **Forward Tilt Option (Standard on Electrically Operated Chairs):** Benefit from the added forward tilt feature for a more ergonomic experience.
- **Electric Option Available:** For those seeking automation and ease of use, the Configura Advance is also available in an electric version.

SPECIFICATIONS

Specification	Details
Colour	Charcoal
Fabric	Vinyl fabric with VP on all contact areas
Seat Width Options	350 mm (14 in) – 550 mm (22 in)
Seat Depth Options	400 mm (16 in) – 550 mm (22 in)
Seat Height Options	330 mm (13 in) – 530 mm (21 in)
Tilt-in-Space	37 degrees
Backrest Recline	130 degrees
Forward Tilt	10 degrees
Weight Capacity	160 kg (25 st)
Mechanism	Manual tilt-in-space

VIDEO

[embed]https://www.youtube.com/watch?v=msXT6VjWgoU[/embed]

BROCHURE

Configura Advance Brochure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au