

CONFIGURA COMFORT & ADVANCE ELATA E600C ALTERNATING AIR CUSHION

SKU: E600C

PRODUCT DESCRIPTION

Purpose & Benefits

The **Elata E600C Alternating Air Cushion** is a dynamic pressure relief solution designed for individuals who have developed pressure injuries. Using alternating inflation and deflation cycles, it redistributes pressure to aid in the prevention and management of skin breakdown.

This cushion is particularly suited to users with pressure injuries in the seating area that have not responded to foam cushions, or those at very high risk due to immobility or the presence of a persistent ulcer. It provides enhanced comfort and accommodates users weighing up to 160 kg.

Features

- Alternating cell inflation and deflation for cyclic pressure redistribution
- Provides comfort and assists with existing pressure injury management
- Compatible with Configura Comfort and Configura Advance chairs

Specifications

- **Includes:** Cushion and pump
- **Maximum User Weight (SWL):** 160 kg

- **Compatibility:** Configura Comfort and Configura Advance chairs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au