

ERGONOMIC FOREARM CRUTCHES

SKU: JAN-124-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

These Ergonomic Forearm Crutches are Ergonomically shaped anatomical handgrip forearm crutch designed to distribute pressure more evenly across the hand.

KEY FEATURES

- Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong.
- Ergonomic palm grip offers pressure distribution for maximum comfort.
- Fitted with high quality rubber tips for stability.
- Both hand to floor and palm grip to arm cuff sections are adjustable.
- Maximum User Weight - 160kg

SPECIFICATIONS

- M (Hand Grip Height - 66 - 94cm, Arm Ring Adjustment- 20.5 -27.5cm)
- L (Hand Grip Height - 83 - 111cm, Arm Ring Adjustment- 20.5 -27.5cm)

Buy the Ergonomic Forearm Crutches Australia online at MediMart using Afterpay or NDIS method today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au