

DAYS 102 ROLLATOR

SKU: DAYS-102-RED

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Days 102 Series, where lightweight meets strength. Crafted with a high-quality aluminium frame, this design offers the perfect balance between a lightweight structure and enduring build quality.

FEATURES

- **Durable 6″ Wheels:** Navigate smoothly on various terrains.
- **Ergonomic Handgrips:** Designed for comfortable and easy gripping.
- **Adjustable Handle Height:** Customise your grip with 5 different settings at intervals of 25mm.
- **Rest Seat:** Take a moment to sit and relax during your outings.
- **Storage Bag:** Conveniently store your essentials on the go.
- **Secure Braking System:** Featuring loop lockable brakes for enhanced security and convenience. Plus, adjust the brake tension effortlessly to fit your needs.
- **Foldable Backrest:** Ensures comfort and easy storage.
- **Easy Fold Mechanism:** No hassle storage and transport. Simply pull the seat strap and fold it.
- **Lightweight:** Weighing only 6kg, it’s easy to manoeuvre and transport.
- **Strong Weight Limit:** Supports a safe working load (SWL) of up to 125kg.

SPECIFICATIONS

Rollator Type	4 Wheel
Frame Material	Aluminium
Seat Height Front	500 mm

Seat Height Rear	500 mm
Seat Width	350 mm
Seat Depth	320 mm
Overall Depth	640 mm
Overall Height	755-875 mm
Gross Weight	6 kg
Max SWL	125 kg
Brand	Days Healthcare

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au