DAYS SHOPPER 8 TRANSIT WHEELCHAIR



SKU: MFI-69397

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

TheShopper 8 Attendant Propelled Wheelchair is ideal if you are a Carer and need to transport a loved one to various destinations as it is lightweight, compact and easy to transport. Constructed out durable aluminium it is ultra lightweight at only 9kg and has an incredible 100kg weight capacity. If safety is important to you then you will love how this wheelchair also comes with a safety lap belt and rear wheel locks.

FEATURES

- Folding Backrest
- Fixed Arm Rests
- Rear Wheel Locks
- Padded Armrests
- Safety Lap Belt
- Attendant Handbrakes
- Height Adjustable, Flip-Up Footplates
- Swing-Away, Leg rests
- Warranty - 12mths Limited warranty.

SPECIFICATIONS

- Weight of Chair - frame only 9kg
- Weight of Chair -with leg rests 105kg
- Overall Seat Width: 432mm
- Seat Width (between arms): 381mm

- Seat Depth: 407mm
- Backrest Height: 432mm
- Leg Rest Height Range: 356 - 483mm
- Armrest Height (seat to top of handles): 254mm
- Wheels (solid castors) Front: 203mm x 25.4mm
- Wheels- Rear: 203mm x 25.4mm
- Maximum User Weight: 100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au