

DAYS SWIFT LITE ATTENDANT PROPELLED WHEELCHAIR

SKU: ANI-ALTS1703-15

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Days Swift Lite Attendant Propelled Wheelchair is ideal for occasional use. The Swift Lite is designed with the carer in mind, the compact lightweight frame makes it easy to transfer in and out of a car. The push button height adjustable armrests make it versatile allowing you to push under a table then quickly adjust for comfort. Attendant propelled in one size.

KEY FEATURES

- Net Weight 9.5kg
- Folding backrest
- Rear wheel locks
- 17” seat width
- Swing-away, removable leg rests
- Height adjustable, flip-up footplates

SPECIFICATIONS

	Days Swift Lite Wheelchair, Attendant Propelled
Overall Height:	895mm
Overall Width:	555mm
Overall Depth:	1050mm
Seat Width:	430mm
Seat Depth:	425mm
Seat Height at Front:	455mm

Seat Height at Rear: 430mm

Buy this Swift Lite Transit Days Wheelchair online with NDIS Payment Method or Afterpay, shop Now!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au