

DAYS WHIRL TRANSIT WHEELCHAIR 18"

SKU: KUN-WHIRL45TR

PRODUCT DESCRIPTION

The DAYS Whirl Wheelchair Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about.

Key Features

- Quick-release wheels (for transport)
- Removable armrests
- Swing away and removable leg rests
- Solid tyres (no maintenance required)

Specifications

- Weight: 16.8kg (complete)
- User Weight: 120kg
- Seat Width: 450mm
- Seat Depth: 430mm
- Backrest Height: 440mm
- Seat Height - Front: 500mm
- Seat Height - Rear: 475mm
- Front Castor Wheel: 8" x 1" PU
- Rear Wheel: 24" x 1 3/8" PU
- Overall Width: 640mm
- Upholstery Nylon (Flame Retardant)
- Colour: Silver

Buy the DAYS Whirl Transit Wheelchair using AfterPay or NDIS method Today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au