

DEUTSCHER FLOOR LINE BED WITH HIGH PROTECTION SYSTEM

SKU: WMR1080HPS

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The High protection system (HPS) has been specially designed to offer maximum comfort and protection to a range of patients suffering from continuous movement, including Huntington’s, Cerebral Palsy, Acquired Brain Injury, Psychosis and Epilepsy. The HPS provides the best possible protection without compromising on the functionality of the Floor line bed operation.

The HPS fits directly onto all Deutscher Floorline beds. The system can be removed easily to convert back to a standard floorline bed.

FEATURES:

•Strong but soft padded side panels and bed ends for user protection.
•Hinged fold down side panels allow greater access to the user.
•The sides can be completely removed, if required for even greater access when utilising patient hoists or when better access to the user is required.
•The complete High Protection System can be removed when not required, converting the bed back to a standard floor level bed.

Please View 'Specifications' Tab above for more information.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au