

DRIVE AUTOMATIC FOLDING SCOOTER

SKU: S3051-red

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The latest advanced automatic folding travel mobility scooter from the Solax stable is top of the class. Solax Ezy Auto Fold features remote control folding function, more leg room, adjustable tiller height and angle, comfortable seat and armrests, lithium battery and lightweight!

FEATURES

- Automatic folding at the press of a button
- Manual folding option
- Remote control
- Lightweight
- Lithium battery
- Easy Operation
- Adjustable tiller height and angle
- Comfortable seat
- Comes with armrests
- Latest Solax technology
- Speed control
- Automatic braking

SPECIFICATIONS

- **Folding Operation:** Automatic With Remote Control
- **Overall Length:** 930mm

- **Overall Width:** 450mm
- **Tyres:** Solid front 6” Rear 7”
- **Maximum speed:** 6km/hr
- **Speed Control:** Yes
- **Safe Gradient:** 0-12 degrees
- **Weight Capacity:** 125kg
- **Range:** 15km
- **Turning Radius:** 1.4m
- **Weight:** 23.4kg
- **Drive System:** Rear Wheel Drive via Transaxle
- **Seat Height From Ground:** 560mm
- **Battery:** Lightweight Lithium 24v 10ah
- **Motor:** 120w
- **Battery Charger:** 24v 2ah
- **Colours:** Red

Buy the Auto Folding Mobility Scooter Australia online at MediMart using Afterpay or NDIS method today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au