

ARJO EASYTRACK FS - FREE STANDING CEILING HOIST GANTRY

SKU: 9220010

PRODUCT DESCRIPTION

Purpose & Benefits of the EasyTrack FS Free Standing Ceiling Hoist Gantry

The **EasyTrack FS Free Standing Ceiling Hoist Gantry** is a portable and flexible solution for safe and efficient transfers in homes or care facilities. Designed to work seamlessly with ceiling hoists, this gantry system eliminates the need for permanent ceiling installation, providing a convenient alternative for users and caregivers who require lifting assistance.

Key Features

- **Free-Standing Design:** The EasyTrack FS requires no permanent installation, offering a freestanding structure that can be set up in various locations. It's ideal for temporary use or where structural modifications to the ceiling are not possible.
- **Sturdy and Stable:** Engineered for safety, the gantry system is constructed with high-quality materials that ensure stability during transfers, providing peace of mind for both users and caregivers.
- **Portable and Adjustable:** The gantry is easy to assemble, disassemble, and transport, making it a flexible solution that can be moved between rooms or locations as needed.
- **Compatible with Ceiling Hoists:** The EasyTrack FS is designed to work with a variety of ceiling hoist models, ensuring a smooth and safe lifting process for users with limited mobility.
- **No Drilling or Fixing Required:** Unlike traditional ceiling hoists, the EasyTrack FS does not require drilling into walls or ceilings, making it a non-invasive and adaptable option for any environment.
- **Spring loaded height adjustment**
- **Automatic locking mechanism**

Specifications

Height	1729mm – 2305mm
Outside Length	2350mm – 3370mm
Inside Length	1907mm – 2927mm
Total Weight	45kg
SWL	182kg

The **EasyTrack FS Free Standing Ceiling Hoist Gantry** is the perfect solution for those in need of a flexible, portable, and safe hoisting system. Its non-invasive design and compatibility with ceiling hoists make it an excellent choice for enhancing mobility and independence in a variety of care settings.

- Folding sling hooks
- Reduced space for storage/transport
- Electronic microprocessor
- Soft-start and stop motor control
- Unit weight including battery: 6.5kg
- Lifting speed: 4.3cm per second
- SWL: 200kg

Safety Features

- Electrical up and down emergency buttons
- Emergency push button stopping device
- Emergency brake
- Manual emergency lowering device

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au