


ELISPOON SELF STABILISING SPOON

SKU: ELI-ELISPOONXL

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

How is ELISpoon different from regular spoons?

ELISpoon is easy to use. Whatever movement of the wrist is performed, the bowl of the spoon remains horizontal and keeps the food in place. The material used is lightweight and strong. A key characteristic is the spoon’s affordability. There are no electrical components or batteries, making it carry anywhere and to clean. ELISpoon is a product designed and built to ease the life of the users and their love.

What is the benefit of using ELISpoon?

Using ELISpoon, the individual gains self-confidence and recovers a sense of control. Everyone should enjoy and share moments of conviviality at the dining table, which in essence brings people together and blurs differences. That’s what Eli is about. Helping individuals achieve an ability to participate in basic commonplace daily activity.

KEY FEATURES

- ELISpoon is a revolutionary spoon that enables independent eating for people with hand coordination difficulties such as Cerebral Palsy, Spinal Cord Injury, Post-stroke deficits, or Huntington's disease.
- Counterweights and rotation axes keep the spoon bowl level regardless of how your hand or arm twists, bends, or moves, preventing food spillage.
- ELISpoon has no electrical components or batteries, making it easy to use and easy to clean.

Buy this EliSpoon Self Stabilising Spoon online at Medimart as we have NDIS payment method available.

SPECIFICATIONS

Dimensions	ELISpoon XL Stabilising Soup Spoon
Total Length	21cm
Total Width	10.5cm
Handle Length	Approx 12cm
Diameter	1.9cm
Weight	78g
Metal Spoon Head	5.5cm (L) x 3.7cm (W)
Colour	White

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au