

ETAC AKTIV REACHER WITH POWER GRIP AND HOOK

SKU: 80505005

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Etac Aktiv reacher with Power Grip and Hook is specifically designed to provide support for individuals with weak hands and limited lifting ability. It features a forearm support that assists in maintaining stability while using the reacher.

FEATURES

- Similar to other Aktiv reachers, the jaws of this model are turnable, allowing for versatile use.
- Equipped with hooks that enable both pushing and pulling of items
- To enhance durability and improve gripping capabilities, the jaws are coated with hardwearing rubber.
- A magnet is incorporated into the reacher's design, enabling it to pick up small metallic items conveniently.

SPECIFICATIONS

Item Length	Item number	Design	Frame material	Lifting capacity, max (kg)	Material	Weight (g)
60 cm	80505005	Ergonomidesign	Aluminium	2	Handle and jaws: Polycarbonate/ABS. Forearm support and jaw bracket: Polyamide	139

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au