

ETAC FEED ADJUSTABLE SPOON

SKU: 8040330-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The ETAC Feed Adjustable Spoon is specially designed for children to facilitate eating by themselves. Using a power grip, the hand is prevented from sliding thanks to the thumb support. The angle is set with a little screw. The elastic strap is useful for children with poor gripping ability. Etac Feed adjustable spoon is available for both right-handed and left-handed use.

FEATURES

- Lightweight design for ease of use
- Available for left or right-handed users
- Dishwasher safe

SPECIFICATIONS

Weight	40 g
Dimensions	16 cm
Material	White polyamide
Design	Ergonomidesign

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au