

EXOFORM CARPAL TUNNEL WRIST

SKU: OSS517073-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Exform Carpal Tunnel Wrist is a lightweight and low profile wrist brace that provides immobilisation to the joint.

FEATURES

- Lightweight low-profile design
- More room at base of the thumb, palmar crease and 5th MP for full MP flexion
- Designed to provide smooth controlled resistance to extreme range of motion
- Pre-shaped into a neutral position for comfortable limited motion/flexibility
- Provides radial and ulnar styloid relief
- Soft touch web space is cushioned and 40% thinner than other brands
- Interior is 3-dimensionally moulded to enhance comfort

SPECIFICATIONS

Wrist Circumference*	**Palm Circumference****	
Small	13 - 16.5cm	17 - 20cm
Medium	16.5 - 20cm	20 - 23cm
Large	20 - 25cm	23 - 27cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au