HI-LO SCOOP PLATE



SKU: PAT-1540

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Hi-Lo Scoop Plate is ideal for the one-handed user. The vertical wall guides food onto utensils. This light Scoop Plate has a non-slip rubber coated base to prevent sliding.

FEATURES

- Ideal for the one-handed user.
- The vertical wall guides food onto utensils.
- This lightweight plate has a non-slip rubbercoated base to prevent sliding.
- Entry wall is 13mm high and steadily increases upward to 38mm.
- 230mm diameter provides plenty of space for any size meal.
- This microwaveable plate has a 19mm rim for easy handling.
- Top rack dishwasher safe.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au