TERRY CLOTH SOCK AID



SKU: PAT-A7273

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Homecraft Terry Cloth Sock Aid makes sock application easy. The flexible, plastic sock aid has slits in the plastic core allowing it to bend around the heel. The inside is lined with nylon to reduce friction. The outside is lined with terry cloth to hold the sock firmly until it is on the foot. Two pulling tabs with loop handles permit one-handed application. Hand wash only.

FEATURES

- Length 230mm
- Weight 73g

Looking for a Medimart alternative?

Discover the **Medimart Flexible Sock Aid** designed to help users put on socks more easily with reduced bending, twisting and strain for improved independent daily living.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au