

ICARE ACTIVEFLOAT TOPPER OVERLAY

SKU: ICOF-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Icare ActiveFloat Overlay, utilising cutting-edge technology, employs a highly resilient and exceptionally breathable substance constructed to mitigate shear forces and provide pressure alleviation. This top-grade overlay minimises shear in all directions, enhancing the user's comfort. Shear reduction is achieved through the implementation of SlideSkin technology, which employs an additional membrane allowing the cover to correspond with the movements of the user.

FEATURES

- Ventilated Side Fabric - The side wall of this overlay cover utilises Vant fabric, a three-dimensional material made of polypropylene that is highly breathable and resistant to bacteria, thereby significantly improving air circulation
- SlideSkin™ technology reduces shear with a second membrane, enabling the cover to move or slide with the users movements
- Suitable for use on some alternating air mattresses. ActiveX™ Overlays are often used to increase comfort and reduce heat loss when used with alternating air mattresses

SPECIFICATIONS

- **Long Single:** 2030mm x 900mm x 50mm - Min. 30kg | Max. 215kg
- **King Single:** 2030mm x 1070mm x 50mm - Min. 30kg | Max. 215kg
- **Short Double:** 1900mm x 1350mm x 50mm - Min. 30kg | Max. 215kg
- **Queen:** 2030mm x 1520mm x 50mm - Min. 30kg | Max. 215kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au