
ICARE FULL LENGTH FOLD-DOWN BED RAILS - PAIR



SKU: ACFSRP

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Full length fold-down bed rail pair. Your choice of the colour finish allows for a clean finish in all environments. Easy installation yet an incredibly strong structure means you can feel safe and secure without sacrificing space, aesthetics or comfort. **For ICARE Beds ONLY.**

- Compatible with all Icare beds
- Fold down design
- 4 Rails for maximum safety
- Slip on attach system
- Position anywhere along to side as needed
- Charcoal colour
- Maximum user weight: 200kg
- Width: 1385mm
- Height above mattress platform: 495mm
- Space between rails: 80mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au