

ICARE BEDS - IC333 HOMECARE BED

SKU: IC333

PRODUCT DESCRIPTION

ICARE IC333 HOMECARE BED BENEFITS

The ICare IC333 is an electrically operated bed that has been designed for both home care and healthcare bed facilities. This bed has adjustable back rest, knee break, is height adjustable and also offers the user trendelenburg tilt. The IC333 also uses a flexible slat systems that support the mattress and comes with an upholstered surround. Available in Onyx or Stone fabric.

Companion Beds and head boards also available.

Key Features

- Maximum load rated 200kg
- Upholstered head board to match IC333
- Back rest and knee rest
- Trendelenburg tilt 15 degrees
- Heavy-duty locking castors
- Choice of fabrics for upholstered
- Dual lift motors totalling 12,000 Newtons lifting power

Size Options

- IC333LS
  Long Single
  2030 x 900
- IC333LS
  King Single

2030 x 1070

- IC333KS
  Long Double
  2030 x 1350
- IC333D
  Queen
  2030 x 1520
- IC333Q
  Split Queen
  2030 x 750 (per side)
- IC333SQ
  Split King (partner set-up)
  2030 x 900 (per side)

Colour & Fabric Options

View the full range of options here - ICare-Upholstery-Colour-Choices.pdf

Accessories

- High side rail
- Low side rail
- Bed sticks
- Full length fold down bed rail
- Over bed table
- Over bed pole

Specifications

- Safe Working Load:
  - 200kg for Long Single & King Single (Includes any Icare mattress)
  - 300kg for Double & Queen (Includes any Icare mattress)
- Height Range:
  - Minimum Height: 22cm (floor to top of mattress platform)
  - Maximum Height: 66cm
- Lift Angles:
  - Head Lift Angle: 62°
  - Knee Break Angle: 42°
- Motors & Actuators:
  - 1x 6000N / 24W Head Adjust Actuator
  - 1x 6000N / 24W Knee Adjust Actuator
  - 2x 6000N / 24W Height Adjust Actuators (2x 8000N for Double & Queen)
  - No-Pull Motors for added safety

- Hand Control Unit:
  - Backlit for Night Use
  - Wired Stretch Cord
  - Hook Attachment for Side of Bed
  - Wipeable for Hygiene
- Castors:
  - 8 Heavy-Duty Swivel Castors
  - One Lockable Castor on Each Corner

Colour & Fabric Options

- Premium Fabric: Bison (Vinyl)
- Premium Fabric: Mist
- Premium Fabric: Noir
- Premium Fabric: Royal (Discontinued)
- Premium Fabric: Sand (Vinyl)
- Premium Fabric: Slate (Vinyl)
- Onyx
- Silver
- Smoke
- Stone
- Vinyl Dove
- Vinyl Ebony

Buy the ICare Bed IC333 at Medimart Australia, as we have NDIS Payment method available.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au