

INVACARE BIRDIE EVO LIFTER

SKU: 1650165-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The new Birdie Evo has been designed to offer a high feature-set of options that have been blended with new technologies to provide maximum comfort and security when lifting or transferring a patient to or from a bed, chair or even the floor.

Based on the renowned Birdie, the new range has been customer insight led, encompassing everything that was valued on the Birdie whilst bringing in new features to further enhance the client and carer experience.

The Birdie Evo encompasses an ergonomic lift that allows the client to be easily rotated 360°, making positioning much easier for the carer and reducing anxiety for the client.

The high boom ensures there is space for the client, even when the hoist is in its highest position, and the optimised space in front of the actuator reduces the risk of injury to the client’s knees.

FEATURES

- Securely attach slings with one hand thanks to a new user-friendly hook design
- Exchange spreader bar quicker, more safely and easier than ever before with the new quick release mechanism. No tools required
- Integrated dampner greatly reduces rocking movement
- Ergonomic push handles

SPECIFICATIONS

- BIRDIE EVO Lifter – 180 kg Capacity
  The BIRDIE EVO Lifter is designed with an electric leg lifter and is equipped with the CBJC system (Control Box with JumboCare Electric).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au