

INVACARE REA AZALEA TRANSIT WHEELCHAIR

SKU: AZALEA-TR-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Invacare Rea Azalea Transit Wheelchair is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes Azalea the ideal choice. The new backrest, Flex 3, is multi-adjustable in height, width, and shape for different body forms. A self-propelled version is also available.

FEATURES

- The Invacare Rea Azalea Transit Wheelchair features a Flex 3 tension adjustable backrest and backrest plates, providing superior customisation to ensure comfort and proper posture.
- The 25° seat tilt and 30° backrest recline allow users to find the most comfortable sitting position, reducing strain and optimising long-term usage.
- Unique to this model is the biomechanical positioning hinge on the backrest, promoting natural movement and increased mobility for the user.
- The wheelchair is available in 3 seating width options, offering a personalised fit to accommodate various body types.

SPECIFICATIONS

- Width: seat width + 23cm
- Length: 102cm
- Product Weight: 34kg
- Warranty: 2 years
- Max User Weight: 135kg
- Seat Width: 40.6cm (16"); 45.72cm (18"); 50.8cm (20")
- Seat Depth: 43cm - 50cm
- Seat Height: 45cm

- Tilt Range: 25°

COMPLIANCE

ARTG 99518
AS|NZS 3695
Crash Tested: ISO 7176-19

WARRANTY

Frame: 2 years | Components: 2 years

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au