

INVACARE SELF HELP POLE - CS8 BED COMPATIBLE

SKU: CS8-SP002

PRODUCT DESCRIPTION

Purpose and Benefits

The Invacare Self Help Pole for the CS8 Bed is designed to assist patients in repositioning themselves while in bed. By providing a stable overhead support handle, it enables users to shift, lift, or adjust their position independently, enhancing comfort and supporting mobility during recovery or long-term care.

This accessory promotes user autonomy, reduces dependence on carers for minor adjustments, and is ideal for hospitals, aged care facilities, and home care environments.

- Supports safe, independent repositioning in bed
- Encourages mobility and active participation in care
- Helps reduce manual handling strain for carers
- Suitable for a wide range of care environments

Features

- Specifically designed for compatibility with Invacare CS8 beds
- Mounts securely to either side at the head-end of the bed

- Durable powder-coated steel frame
- Easy to install and remove as required
- Contoured grip handle for secure hand positioning
- Complies with Australian safety and quality standards

Specifications

- Model: Self Help Pole – CS8 Bed
- Mounting Position: Left or right side at head-end of bed
- Compatibility: Invacare CS8 Single and King Single beds
- Construction: Steel frame with powder-coated finish
- Handle: Angled with textured grip
- Weight Capacity: Up to 75 kg
- Warranty: 1 year

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au