JUDITTA TILT IN SPACE WHEELCHAIR



SKU: 1008-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Juditta is a superb wheelchair and ergonomic seating system that provides stable, comfortable positioning. Especially recommended for individuals with different levels of disability; it even adapts to the most complex cases. Ideal for clients who need an extended or permanent stay in wheelchair or bed. The ergonomic design supports and holds body shape to help prevent contractures and re-distribute pressure. Weight is well distributed between the seat and backrest which helps avoid pressure peaks on the ischial area. This level of comfort helps increase sitting position tolerance and provides correct alignment of the trunk, without creating pressure on the spine. The ergonomic height-adjustable backrest has a unique structure which provides correct alignment to the trunk and a high level of stability without the need of side supports. The Tilt & Recline feature makes it possible to obtain the ideal support position. When you change the tilt of the leg rest and back rest, hip and knee joints will not change. It also positively affects blood pressure by improving oxygenation in the skin and muscle tissue. Tilt & Recline functions are operated by two independent gas springs that absorb body weight and make it easy to be handled by the caregiver. Unique semi-elastic cords are highly breathable, allowing air to flow for an optimal micro-climate. This design adapts to any type of body shape and weight and provides considerable reduction in pressure peaks. Even after one hour from positioning, the seating cords continue to maintain their shape and pressure loading decreases. The special seat and back cover is made from fire-retardant material which is highly breathable, tough (maximum resistance to rubbing), yet breathable and anti-slip (an exclusive Ormesa formula). Seating covers come in long-lasting and highly resistant colours; they are resistant to sweat, light and can be dry cleaned or washed in washing machine. The multi adjustable headrest provides stability, without obstructing the lateral movement of the head. The headrest follows the shape of the backrest, thus assuring the proper alignment between head and trunk. Adjustable in height and

depth. All adjustments improve the visual orientation. Made from anti-bacterial material with closed cells highly resistant to external chemical agents.

FEATURES

- Models B30 & B60 are crash tested according to ISO 7176/19 and ANSI RESNA WC/19 standards
- Multi adjustable headrest
- Ergonomic backrest
- Easy to transfer
- Hygienic and Breathable
- Comfortable and long lasting seat and backrest

For more information please check out the brochure below.

Click here to view brochure.

SPECIFICATIONS

JUDITTA B12	**SIZE 36**	**SIZE 40**	**SIZE 45**	**SIZE 50**
Seat bottom width	36 cm	40 cm	45 cm	50 cm
Seat bottom depth	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm
Backrest width	47 cm	47 cm	52 cm	57 cm
Overall width	60 cm	60 cm	65 cm	70 cm
Seat bottom inclination-tilt	da 10° a 35°	da 10° a 35°	da 10° a 35°	da 10° a 35°
Backrest inclination	da 10° a 38°	da 10° a 38°	da 10° a 38°	da 10° a 38°
Total weight	37 kg	37 kg	39 kg	40 kg
Max capacity	120 kg	120 kg	130 kg	140 kg
User Height	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm

JUDITTA B30	**SIZE 36**	**SIZE 40**	**SIZE 45**	**SIZE 50**
Seat bottom width	36 cm	40 cm	45 cm	50 cm
Seat bottom depth	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm
Backrest width	47 cm	47 cm	52 cm	57 cm
Overall width	62 cm	62 cm	67 cm	72 cm
Seat bottom inclination-tilt	da 10° a 35°	da 10° a 35°	da 10° a 35°	da 10° a 35°
Backrest inclination	da 10° a 38°	da 10° a 38°	da 10° a 38°	da 10° a 38°
Total weight	38 kg	38 kg	41 kg	42 kg
Max capacity	120 kg	120 kg	130 kg	140 kg
User Height	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm

JUDITTA B60	**SIZE 36**	**SIZE 40**	**SIZE 45**	**SIZE 50**
Seat bottom width	36 cm	40 cm	45 cm	50 cm
Seat bottom depth	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm	40, 43, 46, 49 cm
Backrest width	47 cm	47 cm	52 cm	57 cm
Overall width	69 cm	69 cm	74 cm	79 cm
Seat bottom inclination-tilt	da 10° a 35°	da 10° a 35°	da 10° a 35°	da 10° a 35°
Backrest inclination	da 10° a 38°	da 10° a 38°	da 10° a 38°	da 10° a 38°
Total weight	40 kg	40 kg	42 kg	43 kg
Max capacity	120 kg	120 kg	130 kg	140 kg
User Height	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm	da 140 a 190 cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au