

KCARE DELUXE BOWLS SEAT WALKER

SKU: KA377BW

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The KCare Deluxe Bowls Seat Walker is designed for outdoor use with ultra-wide wheels that can be used on outdoors lawn bowls greens and other outdoor environments. It is height adjustable to suit all users of various size. It has Soft Touch locking brakes and comes with a bag and basket. Also contains an extra wide polyurethane seat. The Deluxe Bowls is the ultimate safety and strength walker to support the user. Its ergonomic design facilitates spacious support for the user.

FEATURES

- Ultimate safety and strength to support the user
- Extra large seat for ultimate comfort
- Adjustable to suit all users of various size
- Ergonomic design to facilitate support for the user
- Easily foldable frame, lightweight aluminium frame
- Soft Touch Braking system
- Extra soft PU tyres as to not damage grass
- Ideal for bowling, golfing, and grass activities
- Privacy pocket under seat

SPECIFICATIONS

Max User Weight	200 kg
Min Handle Height	780 mm
Max Handle Height	920 mm

Overall Width	620 mm
Seat Width	410 mm
Seat Depth	350 mm
Floor To Seat Height	520 mm
Diameter	200 mm
Product Weight	6 kg
Brand	K Care

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au