

KCARE ECONOMY PADDED STEEL SHOWER STOOL HEIGHT ADJUSTABLE

SKU: 12325S

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Economy Shower Stool with Padded Seat** provides a safe and comfortable solution for users who experience difficulty standing while showering. Its padded PU seat offers extra comfort, while the durable steel frame ensures strength and stability. This stool promotes independence and reduces the risk of slips and falls in the bathroom.

FEATURES

- Assists individuals who have difficulty showering independently
- Soft PU padded seat provides enhanced comfort and support
- Height-adjustable legs to suit a wide range of users
- Powder-coated 22mm steel tubing for durability and corrosion resistance
- Compact design ideal for smaller bathrooms
- Lightweight and easy to reposition when needed

SPECIFICATIONS

- **Min Overall Height:** 680 mm
- **Max Overall Height:** 810 mm
- **Seat Width:** 410 mm
- **Seat Depth:** 300 mm
- **Seat Height Range:** 460 – 600 mm

- **Max User Weight:** 110 kg
- **Frame Material:** Powder-coated 22mm steel tubing
- **Seat Material:** PU padded

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au