

# KCARE ADJUSTABLE HEIGHT BED CRADLE

**SKU:** KA545

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

The KCare adjustable height bed cradle is designed to fit under the mattress to form a support to keep bedclothes clear of the feet and legs. The cradle is constructed so that one side is placed underneath the mattress and the other positioned horizontally above the bed holding the weight of the bed covers off the body.

Height adjustable from 480 to 700 mm, and made from powder coated steel, this bed cradle ensures you can adjust to suit both the bed and users' needs.

## FEATURES

- Quick and easy to fit or remove
- Fits under the mattress
- Keeps bedding & clothes clear of feet and legs

## SPECIFICATIONS

| | |
|---|---|
| Overall Width | 510 mm |
| Overall Length | 310 mm |
| Min Overall Height | 480 mm |
| Max Overall Height | 700 mm |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au