

KCARE FOREARM WALKING TUTOR

SKU: KA360AG

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed to tutor the user into an upright gait increasing balance and stability .
Height adjustable forearm gutter armrests .
Depth adjustable hand grips .
Folds easily for storage and transportation .
Forearm cuffs rotate laterally and adjust vertically for user comfort and fit .
Sturdy frame provides stable support .

SPECIFICATIONS

Product Weight [kg]	7
Max User Weight [kg]	100
Max Handle Height [mm]	1350
Min Handle Height [mm]	850
Brand	KCare
Colour	Grey
Style	Forearm Support - Rear Glides

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au