

KCARE HANDMADE OPEN SEAT BLUE VINYL

SKU: KA11HOP-B-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **KCare Handmade Open Seat** in blue vinyl by KCare offers robust, hygienic, and customizable support for commode and shower commode applications. Designed to resist everyday wear, this single-density vinyl seat is waterproof, anti-fungal, anti-bacterial, and flame-retardant—making it easy to clean, durable, and suitable for high-use environments. Available in open, closed, left-side open, and right-side open configurations to suit individual needs.

Features

- Single-density vinyl cover—waterproof, easy to clean
- Anti-fungal and anti-bacterial surface for enhanced hygiene
- Flame-retardant and abrasion-tested for safety and longevity
- Superior tear strength for durability
- Multiple opening styles available: open front, closed, left-side, and right-side

Specifications

Code	Description	Seat Width (mm)	Seat Depth (mm)	Hole Width (mm)	Hole Depth (mm)	Hole Position (mm)	Front Opening (mm)	Seat Thickness (mm)
KA11HCOP-B	Handmade Paed, Blue Vinyl	320	320	130	170	80	70	55
KA11HJOP-B	Handmade Junior, Blue Vinyl	400	435	175	225	100	80	55
KA110PU	Polyurethane Blue	445	435	200	260	90	80	40
KA11HOP-B	Handmade Open, Blue Vinyl	445	435	200	260	90	80	55
KA11HRP-B	Handmade RHS Open, Blue Vinyl	445	435	200	260	90	80	55
KA11HLP-B	Handmade LHS Open, Blue Vinyl	445	435	200	260	90	80	55
KA11HRP50-B	Handmade RHS Open, Blue Vinyl	500	435	200	260	100	80	55
KA11HLP50-B	Handmade LHS Open, Blue Vinyl	500	435	200	260	100	80	55
KA11055-B	Handmade Open, Blue Vinyl	550	460	200	260	100	90	55
KA10131	Pressure Relieving, Dartex	445	435	185	225	100	90	65
KA11050PR	Pressure Relieving, Dartex	500	435	185	225	100	80	65
KA11055PR	Pressure Relieving, Dartex	550	460	185	225	100	80	65
KA11060PR	Pressure Relieving, Dartex	600	470	185	225	100	80	65
KA11065PR	Pressure Relieving, Dartex	650	480	185	225	100	80	65

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au