

KCARE PADDED SHOWER STOOL WITH ARMS AND PADDED SEAT

SKU: KA22210

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed for use in the shower, this KCare Padded Shower Stool has a padded backrest and seat. The padded seat is removable for cleaning. The seat height and angle are adjustable. Arms are provided for support while sitting or transferring. Rubber tips on each leg provide stability. The shower stool has an epoxy powder coated zinc treated mild steel, white frame and water proof polyurethane seat and back rest.

FEATURES

- Provides a rest for the user while showering and is ideal in a smaller shower area
- Adjustable seat height and angle facilitates a comfortable seating position
- The clip on seat allows for easy assembly, cleaning and replacement
- Padded Seat option also available
- This shower stool also has the addition of a fixed padded backrest
- The backrest is permanently fixed to the stool frame
- Zinc treated mild steel frame with a powder coated finish

SPECIFICATIONS

Max Seat Height	640mm
Max User Weight	110kg
Min Seat Height	520mm
Overall Depth	455mm
Overall Width	510mm

Product Weight	4kg
Seat Depth	310mm
Seat Width	360mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au