

KCARE PREMIUM OVER TOILET FRAME - ALUMINIUM

SKU: 12247A45-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **R&R Premium Over Toilet Frame – Aluminium** provides stability, comfort, and safety for individuals who experience difficulty using a standard toilet. Designed to assist larger users, it features a durable aluminium frame with a moulded seat and splash guard. Its lightweight yet extra-strong construction ensures reliable support for daily use, while the adjustable legs allow for customised height to suit individual needs. Ideal for use in homes, hospitals, and aged care facilities, this frame promotes independence and confidence in toileting.

FEATURES

- Ultra-light yet exceptionally strong aluminium frame
- Constructed from 32 mm anti-corrosive aluminium tubing with white powder coating
- Height adjustable legs for user comfort and convenience
- Comfortable moulded seat with splash guard and lid
- Designed to meet Australian standards for quality and safety
- Available in 450 mm, 520 mm, and 600 mm seat width options
- Supports users up to 250 kg

SPECIFICATIONS

450 mm Width Model

- **Product Width:** 590 mm
- **Product Depth:** 490 mm
- **Product Height:** from 630 mm

- **Seat to Floor Height:** 480 – 640 mm
- **Seat Width:** 520 mm
- **Seat Depth:** 410 mm
- **Weight Capacity:** 250 kg
- **DVA RAP Schedule:** BE06
- **DVA RAP Code:** 12247A52

520 mm Width Model

- **Product Width:** 590 mm
- **Product Depth:** 490 mm
- **Product Height:** from 630 mm
- **Seat to Floor Height:** 480 – 640 mm
- **Seat Width:** 520 mm
- **Seat Depth:** 410 mm
- **Weight Capacity:** 250 kg
- **DVA RAP Schedule:** BE06
- **DVA RAP Code:** 12247A52

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au