

KCARE PREMIUM PADDED SHOWER STOOL 520MM

SKU: 12325A52

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Premium Padded Shower Stool (520mm)** is designed to provide comfort, safety, and stability for individuals who experience difficulty standing while showering. Its soft PU padded seat offers superior comfort, while the durable 32mm aluminium frame ensures long-lasting strength and resistance to corrosion. Ideal for everyday use, this stool promotes independence and confidence in the bathroom.

FEATURES

- Assists users who have difficulty showering independently
- Soft PU padded seat provides enhanced comfort and support
- Height-adjustable legs for personalised fit and optimal stability
- Constructed from 32mm powder-coated aluminium tubing for strength and durability
- Lightweight, rust-resistant frame suitable for wet environments
- Available in 450mm, 520mm, and 600mm frame widths

SPECIFICATIONS

- **Min Overall Height:** 670 mm
- **Max Overall Height:** 810 mm
- **Seat Width:** 520 mm
- **Seat Depth:** 300 mm
- **Seat Height Range:** 480 – 640 mm
- **Max User Weight:** 250 kg

- **Frame Material:** Powder-coated 32mm aluminium tubing
- **Seat Material:** PU padded
- **Brand:** KCare

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au