

# KCARE TRANSIT MOBILE SHOWER COMMODE WITH SWINGAWAY FOOTRESTS

**SKU:** KA114S-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

The KCare Transit Mobile Shower Commode with Swing away Footrests is a versatile and durable product designed for individuals with mobility challenges. This commode can be used as a shower chair, toilet chair, or transport chair, providing comfort and convenience in one easy-to-use device.

The commode features swing away and removable footrests that can be adjusted to accommodate various leg lengths. The padded seat and backrest provide extra comfort, while the 3-strap backrest ensures secure positioning. The integrated push-handle allows for easy manoeuvrability, and the swing-up arms allow for unrestricted side transfer.

The commode is designed to be used with K Care patented removable seats and is available in various seat widths. It has four 125 mm total locking castors for added stability and safety.

## FEATURES

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind.
- 3 strap backrest for comfort and easy access.
- Integrated push-handle for manoeuvrability.
- Swing-up arms for unrestricted side transfer.
- Swing away and removable height adjustable footrest.
- 4x 125mm total locking castors.
- Designed to be used with a K Care patented removable seats.
- **Seat, pan, and pan carrier not included - sold separately**

There is also a height adjustable version available, click to view the KCare Height Adjustable Mobile

Shower Commode.

## SPECIFICATIONS

| Seat Width | Overall Depth | Overall Height | Seat Depth | Overall Width | Floor To Seat Height | Armrest Type | Armrest Height | Backrest Height | Footrest Type | Product Weight | Max User Weight |
|---|---|---|---|---|---|---|---|---|---|---|---|
| 320 mm | 750 mm | 1005 mm | 320 mm | 485 mm | 550 mm | Loop | 285 mm | 400 mm | Swing Away | 10 kg | 120 kg |
| 400 mm | 905 mm | 1005 mm | 430 mm | 480 mm | 550 mm | Loop | 285 mm | 400 mm | Swing Away | 10 kg | 200 kg |
| 445 mm | 890 mm | 1005 mm | 430 mm | 520 mm | 550 mm | Loop | 285 mm | 400 mm | Swing Away | 10 kg | 200 kg |
| 500 mm | 930 mm | 1005 mm | 430 mm | 630 mm | 550 mm | Padded Loop | 315 mm | 400 mm | Swing Away | 11 kg | 200 kg |
| 550 mm | 930 mm | 1005 mm | 460 mm | 665 mm | 550 mm | Padded Removable | 315 mm | 405 mm | Swing Away | 12 kg | 300 kg |
| 600 mm | 930 mm | 1005 mm | 470 mm | 715 mm | 550 mm | Padded Removable | 315 mm | 405 mm | Swing Away | 13 kg | 300 kg |
| 650 mm | 930 mm | 1005 mm | 480 mm | 765 mm | 550 mm | Padded Removable | 315 mm | 405 mm | Swing Away | 14 kg | 300 kg |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au