

HOMECRAFT INCURVE PLATE SURROUND

SKU: PAT-AA5662Y

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Incurve Plate Surround is innovatively designed to facilitate single-handed eating on any plate. The guard features three clips with the help of which it can be easily attached to any standard plate.

Attach the guard to the plate before placing the food. Once attached, the guard creates an inward sloping face to facilitate food collection and prevent spillage, making it much easier for the user to eat with a single hand.

FEATURES

- This flexible guard can be fitted to ordinary plates to assist with one-handed eating.
- Three clips attach it to the rim and create a high inward sloping face to help with food collection and prevent spillage.
- Plate size 190 to 254mm external diameter.
- Dishwasher safe. Melamine material dishwasher safe to 180 degrees F (82°C).
- Weight: 45g.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton
SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave,
Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella
SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au