

LUXA HOMECARE BED

SKU: LUXA-P

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Luxa represents the apex of Homecare bed design, offering exceptional comfort and care, all the while preserving an individual's dignity and augmenting the aesthetics of their bedroom.

This bed comes equipped with a multitude of high-quality features. These include a backrest with auto-regression capabilities, a manual vascular support system, and an adjustable height ranging from a lofty 725mm to a mere 100mm. Additionally, it provides the possibility to integrate Bluetooth handsets.

Moreover, it is capable of being upgraded with several optional accessories. Among these are a robust horizontal rail, provisions for the inclusion of self-help poles, lights underneath the bed, and many more.

FEATURES

- Extremely Low Minimum Position
- Vertical Lift
- Auto-Regression Backrest
- Ash or Ocean Padded Headboard with full
- Padded Surround

SPECIFICATIONS

Bed Type	Overall Dimensions (mm)	Sleeping Surface (mm)	SWL (kg)	Low/High Position (mm)	Back Angle (degrees)	Knee Angle (degrees)	Reverse Angle (degrees)	Trendelenburg Angle (degrees)

Bed Type	Overall Dimensions (mm)	Sleeping Surface (mm)	SWL (kg)	Low/High Position (mm)	Back Angle (degrees)	Knee Angle (degrees)	Reverse Angle (degrees)	Trendelenburg Angle (degrees)
Luxa Homecare	2160 x 1000	2000 x 900	250	150/725	62	35	15	15
Luxa King Single	2160 x 1170	2000 x 1070	250	150/725	62	35	15	15
Luxa Double	2160 x 1470	2000 x 1370	250	170/725	62	35	15	15
Luxa Queen	2160 x 1600	2000 x 1500	250	170/725	62	35	15	15

BROCHURE

Click here to view the brochure Luxa-by-Caremed-Alrick

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au