

MERITS P312A DUALER POWER CHAIR

SKU: P31231WRMA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Easily turn around in small places with the Dualer Power Chair with an elevating seat. The turning radius is 21” and easy to control by using the swing-away joystick. This versatile power chair reaches the top speed of 8 kmph, has a weight capacity of 136 kilos, and has remarkable handling.

FEATURES

- Easy to convert from RWD to FWD, 2-in-1 power chair, ideal chair for indoor as well as outdoor use.
- Unique compact size for maximum manoeuvrability in narrow space and walkway.
- Height adjustable flip-up armrests can provide seat width adjusted.
- Start with RWD and move to FWD by rotating the chair 180°, to enjoy direct-drive, direct-response feeling.
- Turning radius is only 53cm/ 21”.
- Low seat height of 18” provides excellent stability.
- Power seat elevator available.

SPECIFICATIONS

Driving Range	28 km / 17.5 mil
Max Speed	8 km/h , 5 mph

Battery Type	12V / U1 x 2pcs
Motor	DC 24V, 1240W (Rating) / 950W (Max)
Max Climbing Angle	10°
Turning Radius	53 cm / 21"
Ground Clearance	5 cm / 2"
Controller	PG VS2 60
Charger	5A off-board

Weight Capacity	136 kg / 300 lb
Total Weight	85 kg / 187 lb
Front wheel	8" (200 x 50) PU tire
Drive Wheel	10" (260 x 85) foam filled tire
Rear wheel	10" (260 x 85) Foam filled tire

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au