

MOLIFT SMART 150 PORTABLE PATIENT LIFTER

SKU: M10200

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Awarded for Design Excellence in Europe, the Molift Smart 150 is a lightweight, portable patient lifter designed to be easily wheeled and transported. It is ideally suited for use in home care, for travelling, transporting in a car or to be easily stored.

With a lifting range of 27 – 168 cm (10.6" – 66.1") the Smart 150 provides an excellent maximum lifting height and at the same time gets low enough to facilitate lifting from the floor.

The 4-point sling suspension is standard and guarantees the patient hoist gives a comfortable lifting position for the user with ample space around the head area. This eliminates any risk of the head bumping into the suspension or the user swinging side to side.

With the combination of sling and suspension, the user will slide into a natural, correct seating position, from lying to sitting and vice versa.

The Smart 150 has a unique design with an excellent centre of balance and together with the small footprint makes it easy to manoeuvre even in tight spaces.

FEATURES

- Comfortable and safe transfer
- 4-point sling bar for more spacious lifting position
- Supplied with a 4-point sling bar (suspension). Sling is not included and must be purchased separately.
- Lifts easily from floor

- Small footprint makes it easy to maneouver in narrow spaces
- Low weight
- Folds and unfolds easily without the need for tools
- Easy and flexible to transport when travelling
- Built in software calculates when the lifter requires servicing
- Electric | manual lowering functions
- Emergency stop
- Includes: hand control | battery | charger | sling bar
- Molift Hard Travel Case For Smart Lifter 150 available
- Molift Soft Travel Bag For Smart Lifter Available

SPECIFICATIONS

Molift Smart 150	
Item number	M10200
Battery	NiMH 14.4 V-2.6 Ah
Battery charge time (hours)	3
Hoisting height max (mm)	1665
Hoisting height min (mm)	210
Hoisting interval (mm)	1410
Hoisting speed 75 kg (mm/sec)	60
Leg height (mm)	110
Length (mm)	1200
Material	Aluminium, steel
Maximum user weight (kg)	150
Rated performance	40 hoists at: 75 kg and 500 mm
Service software	Molift Service Tool
Weight (kg)	26
Width chassis, legs open inner (mm)	925
Width chassis, outer (mm)	660

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au