

NOVIS PROCAIR PRIME ADVANCED MATTRESS REPLACEMENT SYSTEM

SKU: APMPC-R03A-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

A highly featured option for those in the high-risk category for pressure injuries, the ProCair® Prime Advanced can be used as a fully alternating pressure mattress or easily converted to a hybrid foam mattress by adjusting the position of the built-in foam layer. It easily and quickly converts back into a dynamic mode when complete pressure off-loading is required. Easy to set up, the ProCair® Prime automatically adjusts the pressure to suit the weight of the patient and does not need a separate pump unit to hang from the bed - it is all built-in! Equipped with additional safety features including back up battery, Auto Fowler Boost, real-time auto weight sensing and cell in cell design, the ProCair Prime Advanced offers additional pressure injury prevention with the same versatility as the ProCair Prime.

The ProCair® Prime Advanced Mattress Replacement System is equipped with an easy-to-use handset that can be operated by the user or the carer. Concerns over correct pressure level settings are eased with AutoCair, a real-time automatic weight-sensing technology that automatically adjusts the pressure settings to keep the client within a safe pressure range for effective therapy, delivering fail-safe patient comfort and support. The inbuilt service and battery life indicators make for worry-free maintenance. Choose from a Single or King Single size from the dropdown menu when making an enquiry.

FEATURES

- **Unique Built-In Pump:** Equipped with a one-of-a-kind built-in pump, which simplifies operation and minimises trip hazards whilst safeguarding against damage that can occur when pumps are accidentally bumped off the bed.
- **Backup Battery:** Providing up to 5 hours of emergency backup power in case of power failure or accidental disconnection.

- **User-Friendly Handset:** Designed for ease of use and can be operated by both the user and the caregiver. Enhances convenience in adjusting settings and modes.
- **"Bottoming Out" Prevention:** Incorporates a lower "safety cell" that remains fully inflated, ensuring greater surface stability and patient safety.
- **HeelCair Technology:** Enhanced protection and therapy for vulnerable heels. Includes an independent heel zone with narrower cells that can be disconnected as needed.
- **AutoCair Technology:** Automatically adjusts pressure settings using continuous weight-sensing technology, eliminating concerns about incorrect pressure levels and ensuring fail-safe patient comfort and support.
- **Real-time Pressure Adjustment:** Automatically adjusts pressure settings in real-time to accommodate weight changes, ensuring effective therapy.
- **Adaptability:** It can function as a fully alternating pressure mattress (dynamic mode). Easily converted to hybrid mode by positioning the foam layer above the air cell layer. Provides flexibility in care and treatment.

SPECIFICATIONS

Single

- Length: 2000 mm
- Width: 900 mm
- Height: 200 mm

King Single

- Length: 2000 mm
- Width: 1050 mm
- Height: 200 mm

Technical Specifications

- Patient Risk Category: High Risk
- Mattress Recommendations: Designed for high-end, home, and aged care needs. Features a unique discreetly integrated pump. Easy to set up with automatic weight sensing. Can be easily converted to a foam hybrid mattress.
- Clinical Features: 2 cell cycle therapy. Full Mattress Replacement. Automatic weight sensing in cycle time. Carer/patient handset controller.
- Quick Twist CPR Valve
- Top Cover: Latex Free

Safety

- Recommended user weight capacity: 30 - 200 kg

Warranty

- 3 years

Compliance

- ARTG Number: 289458
- Australian Standards: IEC 60601-1, IEC60601-1-2, IEC60601-1-11, BS7175 Crib 5

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au