

OSCAR EASY GLIDE DINING CHAIR

SKU: MAEA-DT

PRODUCT DESCRIPTION

OSCAR EASY GLIDE

Dining is made easy for everyone. Lightweight with a simple weight-activated gliding action, this Oscar Easy Glide dining chair increases independence and reduces injuries. Whether you’re in a commercial or domestic setting you can’t go past the Easy Glide.

FEATURES

- Allows the user or carer to effortlessly roll into the table.
- 50mm Rear self-locking wheels (release when seated locks when you stand for stability).
- Light and Rigid Aluminium Frame.
- 38mm front fixed wheels.
- Easy clean wood grain finish.
- 120kg weight capacity.
- High resilience foam.
- Available with matching arm pads.

SPECIFICATIONS

- Seat Height: 480mm
- Seat Width: 500mm
- Seat Depth: 450mm
- Chair Height: 830mm

- Back Height: 370mm
- Chair Width: 570mm
- Arm Height: 690mm

WARRANTY

2 Years Warranty: Mechanism, fabric and electronics

5 Years Warranty: Timber Frame and Foam.

Have a look at other utility chairs.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au