

OSSUR EXOFORM WRIST BRACE

SKU: OSS507087-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Exoform Wrist is designed to provide superior fit and comfort with contoured padding and to provide great support and immobilisation with rigid outer plastic contouring the wrist

FEATURES

- Ergonomically contoured padding on the inside of the brace provides the best fit and most comfort
- Rigid plastic supports contour to the hand and wrist at the top and bottom leading to greater support and immobilisation
- Adjustable aluminium palmar stay allows for fine-tuning when needed

INDICATIONS

- Mild to moderate sprains and strains - Tendonitis - Post-cast healing - Other selected soft tissue injuries

SIZING

Small 16 - 18cm

Medium 18 - 20cm

Large 20 - 23cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au