

OTTOBOCK MOTUS XXL WHEELCHAIR

SKU: 480F61=2_AA05_V01-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Ottobock Motus XXL Wheelchair stands out in the realm of adaptive wheelchairs. It's not just about mobility; it's about delivering a reliable and efficient experience, built around the user's specific needs. With the ability to customise according to preferences and requirements, the Motus XXL promises optimal comfort and enhanced mobility for its users.

FEATURES

- **Robust Folding Mechanism:** Its unique internal cross-brace folding mechanism is both reliable and durable. This ensures ease of storage and transport, making it perfect for those on the go.
- **Enhanced Weight Support:** With the added option of a double cross brace folding mechanism, the wheelchair can comfortably support a user weight of up to 140kg.
- **Constructed with Excellence:** Made from lightweight aluminium, the chair promises sturdiness without compromising on manoeuvrability.
- **Protection & Comfort:** Integrated side panels protect the user's clothing from any contact with the wheels. Moreover, the breathable backrest upholstery ensures the user remains comfortable, even during prolonged periods of sitting.
- **Ease of Transfer:** Swing-away footplates, available as an option, enable effortless transfers and improved access in front of the chair. It also supports foot propulsion, adding an extra layer of mobility.
- **Travel-friendly Design:** The optional half-folding seat back feature allows the chair to fold forward, reducing its overall height. This makes the Motus XXL an excellent choice for travellers, ensuring it can be stored compactly.

SPECIFICATIONS

Frame material	Aluminium
Weight from	approx. 13kg (CV)
Maximum load capacity	125kg (up to 48cm seat width) 140kg (from 50.5cm seat width)
Centre of gravity adjustment	48 positions through the rear wheel adaptor
Wheel locks/brakes	Knee-level wheel lock Drum brake
Seat width	45.5cm (17.9”); 50.5cm (19.8”)
Seat depth	36 - 54cm
Overall length	83 - 111cm
Overall width	52 - 72cm with standard rear wheels and 0° camber 55.5 - 75.5cm with rear wheels with drum brake, 0° camber
Back	Straight, height adjustable / 30° angle adjustable
Back height	30 - 50cm
Lower leg length	34 - 52cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au