

DAYS OVERBED TABLE

SKU: JAN-758KD

PRODUCT DESCRIPTION

PURPOSE BENEFITS

The Days Overbed Table has an improved design featuring an epoxy coated frame with a vinyl top. It also has recessed edges, which enables the table to be cleaned easily. The 50mm twin wheel castors allow the table to be easily moved away from the bed when not in use.

FEATURES

- Height adjustable with turn handle knob
- Epoxy coated frame
- Easy clean vinyl table top
- Recessed edge
- 2 lockable castors

SPECIFICATIONS

- Length: 800mm
- Overall Width: 400mm
- Tabletop Size: 80mm x 40mm
- Minimum Height: 710mm
- Maximum Height: 1050mm
- Height of Castor Wheels: 50mm
- Max Load: 15kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au