P5500 PULL-APART BED



SKU: P5500-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Caremed Alrick P5500 Pull-Apart Bed is a highly functional and versatile medical bed designed to cater to the needs of patients and healthcare providers in various medical settings. This bed is specially crafted to offer ease of use, comfort, and practicality, particularly in environments where space optimization and patient mobility are key considerations.

FEATURES

- **Pull-Apart Functionality:** The unique pull-apart feature allows the bed to be easily separated into two parts, facilitating efficient transport and storage. This is particularly useful in settings where space is limited or the bed needs to be moved frequently.
- **Adjustable Height:** Equipped with an electric height adjustment mechanism, the bed can be raised or lowered to suit the specific needs of patients and caregivers, enhancing comfort and facilitating easier patient transfers.
- **Robust Construction:** Built with durable materials, the P5500 is designed to withstand regular use in a healthcare environment, ensuring longevity and reliability.
- **User-Friendly Controls:** The bed comes with intuitive controls that make adjustments simple and efficient, reducing the workload for caregivers and improving the overall patient experience.
- **Safety Features:** Includes essential safety elements like side rails to prevent patient falls and lockable casters to secure the bed in place during use.
- **Comfortable Sleeping Surface:** The bed is designed with a comfortable mattress support platform, suitable for various types of medical mattresses, ensuring patient comfort during recovery or treatment.

The Caremed Alrick P5500 Pull-Apart Bed is an excellent choice for medical facilities looking for a bed that offers flexibility, comfort, and ease of use. Its innovative design and robust construction make it a valuable asset in providing high-quality care and enhancing the efficiency of healthcare operations.

SPECIFICATIONS

- **Size Specifications:** King Single
- Overall Width - 1070mm
- Overall Length - 2190mm
- Sleeping Surface - 1050 x 2000mm
- Height - 295-740mm
- Safe Working Load - 200kg

- **Size Specifications:** Single
- Overall Width - 930mm
- Overall Length - 2190mm
- Sleeping Surface - 900 x 2000mm
- Height - 295-740mm
- Safe Working Load - 200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au