

REDGUM BARIATRIC PRIMO 8" WALKER

SKU: RG4213GY

PRODUCT DESCRIPTION

Purpose & Benefits of the Redgum Bariatric Primo 8" Walker

The Redgum Bariatric Primo 8" Walker is a heavy-duty mobility aid designed to provide superior support and stability for bariatric users. With its robust construction and ergonomic design, this walker is ideal for individuals who require extra strength and durability in their mobility aids. The height-adjustable seat and handles ensure a comfortable fit for each user, while the folding design makes it easy to transport and store. Suitable for both indoor and outdoor use, the Primo 8" Walker supports users in maintaining their independence and mobility.

Key Features

- **Heavy-Duty Construction**: Built with a strong, reinforced frame, the walker can support users up to 200 kg, providing reliable stability and durability.
- **Height-Adjustable Seat and Handles**: Customise the seat and handle heights to ensure a comfortable and ergonomic fit for various user needs.
- **Large 8" Wheels**: Equipped with 8-inch wheels that offer smooth navigation over different surfaces, ensuring ease of use indoors and outdoors.
- **Comfortable Padded Seat and Backrest**: Features a padded seat and backrest, providing a comfortable place to rest during extended use.
- **Folding Design for Easy Transport**: The walker folds easily for convenient storage and transport, making it ideal for travel or saving space at home.

Specifications

- **Weight Capacity**: 200 kg
- **Seat Height Range**: Adjustable from 55 cm to 62 cm

- **Handle Height Range**: Adjustable from 82 cm to 102 cm
- **Wheel Size**: 8-inch diameter
- **Frame Material**: Reinforced aluminium
- **Folded Dimensions**: Compact for easy storage and transport
- **Cleaning**: Wipe clean with a damp cloth or mild detergent

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au