

REDGUM LIGHTWEIGHT TRI WALKER

SKU: RG4402-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Step into a world of enhanced mobility and daily convenience with the Redgum Lightweight Tri Walker. This modern and sleek design is not only built for stability and comfort but also offers practicality for both daytime errands and night-time strolls.

FEATURES

- **Sturdy Frame:** Crafted with epoxy powder-coated aluminium, ensuring durability and longevity.
- **Contemporary Design:** Modern aesthetics that blend seamlessly with function.
- **Optimal Mobility:** Large wheels coupled with a swivel front wheel make navigating smoother than ever.
- **Customisable Comfort:** Enjoy height-adjustable handles and ergonomic rubber hand grips tailored to your needs.
- **Safety First:** With easy-to-engage/release brake levers, locking handbrakes, and reflective sides/front, you're ensured a secure experience.
- **Additional Perks:** Handy holder for your walking stick ensures you have everything you need at arm's length.
- **Hassle-free Assembly:** Set it up in no time, sans tools!
- **Compact Storage:** Designed to fold effortlessly, making storage and transport a breeze.

SPECIFICATIONS

- Overall Dimensions: Length: 580mm, Width: 675mm, Height: 830mm (lowest setting)
- Handle Heights: 830 - 930mm: 5 settings (floor to top of handle)

- Hole 1: 830mm
- Hole 2: 855mm
- Hole 3: 880mm
- Hole 4: 905mm
- Hole 5: 930mm

- Handle Width: 480mm (overall)
- Wheels: 190mm Diameter
- Product Weight: 4.6kg

- **Please note:** Accessory bag is not included.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au