

# REDGUM SPACE SAVER BOWL AND LID ACCESSORY

**SKU:** RG5BLH

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

A portable or permanent use Toilet Seat Raiser Accessory from the Redgum Space Saver range.

## FEATURES

- To suit the Redgum Space Saver Toilet Seat Raiser (sold separately).
- Converts the Space Saver Toilet Seat Raiser into a Space Saver Bedside Commode.
- Easy to install and remove.
- Easy to clean.
- Lightweight.
- Holds up to 3 litres.

## SPECIFICATIONS

- Overall Dimensions: 310mm (L) x 280mm (W) x 140mm (H)
- Product Weight:
- Bowl: 250g
- Lid: 70g
- Max. Load: 3 Litres

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au