

REGENT COMFORT SLING

SKU: SUPSFB

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Excellent back and neck support.
- Commode opening with divided legs.
- For patients who need full body and neck support.
- Ideal for amputees.
- Available in solid polyester and mesh variants. Polyester mesh sling is ideal for bathing patient.
- Positioning straps.
- Colour coded lifting straps.
- Made to Australian standards.
- N.A.T.A.* Tested safe working load is 250kg.

For sizing see dimensions (cm) in the diagram below:

SIZE

Small A: 1100 B: 650 C: 700

Medium A: 1150 B: 700 C: 750

Large A: 1200 B: 800 C: 800

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au