

REGENT GENERAL PURPOSE SLING

SKU: SUPSGP

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Divided legs can be crossed in the centre, uncrossed or positioned under thighs
- Suitable for patients with good head control
- Available in polyester solid and mesh variants.
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards
- N.A.T.A.* Tested safe working load is 250kg

For sizing see dimensions (cm) in the diagram below:

SIZE

- X-Small: A: 750 B: 900 C: 500 D: 300
- Small A: 800 B: 1150 C: 550 D: 400
- Medium A: 900 B: 1300 C: 600 D: 450
- Large A: 950 B: 1350 C: 650 D: 500

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au