

REGENT TOILET / STANDING SLING

SKU: SUPSTS

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Soft padding throughout sling
- Allows for safe and comfortable chair to toilet transfer
- Suitable for patients who generally have good trunk and hip control
- Allows patient to be easily brought to a standing position
- Elasticated neoprene waist strap
- Easy access for clothing removal during toileting
- Available in solid polyester only
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards
- N.A.T.A.* Tested safe working load is 250kg

SIZE

- Small/Medium A: 965 B: 540 C: 380
- Large/X-Large A: 1080 B: 620 C: 380

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au